



PRODUCTS CATALOGUE





About us.....	Page 2
Our products.....	Page 3
Masks.....	Page 4
Adult Disposable Mask, Kids Disposable Mask.....	Page 5
Medical Mask, Surgical Mask.....	Page 6
FFP2/ KN95 Mask, Face Shield.....	Page 7
N95 Cup , N95 Foldable.....	Page 8
Apparel.....	Page 9
General Isolation Gown.....	Page 10
Isolation Gown AAMI Level 2, Isolation Gown AAMI Level 3.....	Page 11
Sterile Surgical Gown AAMI Level 3, Shoe Cover.....	Page 12
Bouffant Caps, Mob Cap.....	Page 13
Gloves.....	Page 14
Nitrile Gloves, Synthetic Vinyl Gloves.....	Page 15
Syring and Needle.....	Page 16

ABOUT US

ALIT Medical is an entity of ALIT Group Limited, an international trading company.

ALIT Medical specializes in sourcing, managing and executing the production and distribution of Personal Protective Equipment across the globe.

We partner with the world's leading manufacturers and have a dedicated quality control team on the ground at all locations, ensuring quality, competitive pricing and optimal service.

With over a decade of experience on the market, our partners have a proven track record of excellence in consistently providing a steady supply of top-quality products. Our customers vary from private and state healthcare institutions to multinational retail groups in North and South America, Europe and the Middle East.



ALITMEDICAL

OUR PRODUCTS

FACE MASK & SHIELD

Our wide range of face coverings varies from the highest quality disposable masks to medical and surgical masks.



APPAREL

Gowns, caps and shoe covers are all made from the best materials for maximum comfort, varying from Level 1 to Surgical Standards.



GLOVES

Our primary focus is on supply of two types of gloves – Synthetic Vinyl and Nitrile.

We source our gloves from the top manufacturers in Vietnam, Malaysia, Thailand and China.



ALITMEDICAL

FACE MASK & SHIELD



ALITMEDICAL

FACE MASK & SHIELD



ADULT DISPOSABLE MASK

3 layers of non-woven fabric including a layer of melt-blown filter to provide an effective protection.

This is a general use face mask and is not a medical device.

Available in blue, pink and white.

Size: 17.5 x 9.5cm

Applications:

Healthcare, general housekeeping, laboratory, industrial purpose / PPE

Certification:

Meets GB/T 32610-2016
FDA approved (product code LYU, Class 1)
CE approved (product code DY95)
Standard EN 149:2001
+A1:2009



KIDS DISPOSABLE MASK

3 layers of non-woven fabric including a layer of melt-blown filter to provide an effective protection.

This is a general use face mask and is not a medical device.

These masks come in blue, pink and white.
Additional prints are available.

Size: 14 x 9cm

Certification:

Meets GB/T 32610-2016
FDA approved (product code LYU, Class 1)

FACE MASK & SHIELD



MEDICAL FACE MASK

Medical grade masks have at least a 95% droplet filtration system and are comfortable to breath in.

Mask consists of 3 layers: 2 of non-woven spunbonded and a melt-blown.

Medical masks are rectangular in shape.

Each layer has fine to very fine fibers – all latex free.

Size: 17.5 x 9.5cm

Applications:

Medical purpose / examination, healthcare, general housekeeping, laboratory, industrial purpose / PPE.

Certification:

EN 14683 standard
95% bacterial filtration efficiency
FDA approved (product code QKR, Class 1)
CE approved (medical mask)
ASTM Level 2



SURGICAL FACE MASK

A surgical mask is intended to be worn by healthcare professionals during various procedures. It is designed to protect patients and treating personnel by catching bacteria shed in liquid droplets and aerosols from the wearer's mouth and nose.

Surgical masks vary in quality and levels of protection.

Despite their name, not all surgical masks are suitable to be used during surgery.

Mask contains 3 layers – two of non-woven spunbonded and a melt-blown.

Size: 17.5 x 9.5cm

Applications:

Medical purpose / examination / surgery.

Certification:

EN 14683 standard offering 98% bacterial filtration efficiency
FDA approved (product code MSH class 2)
CE approved (surgical mask)

FACE MASK & SHIELD



FFP2 / KN95 FACE MASK

KN95 respirator mask designed to capture 95% of 0.3-micron particles.

These respirator face coverings are made of 5ply filter layers for extra protection. This mask can last all day offering a comfortable fit.

Protection from a variety of airborne particles, including those found outdoors, droplets and allergens.

Size: 10.15 x 15 x 25cm

Applications:

Healthcare, general housekeeping, laboratory, industrial purpose / PPE

Certification:

China: GB 2626-2006 (type KN95)

FDA approved (product code QKR / MSH, Class 1)

CE approved, standard EN 149:2001 +A1:2009



FACE SHIELD

Face shields are an effective way to protect against airborne droplets and splashes of liquid. They are essential for healthcare and emergency services professionals, as well as other key workers. Face shields are anti-fog PET plastic, which protects against liquids and ensures you can clearly see through it. The face shields are fitted with a foam pad for a comfortable extended use.

Size: 330 x 229mm

Certification:

Meets ANSI/ISEA Z87.1-2015

Meets GB14866-2006

CE approved

FDA approved

FACE MASK & SHIELD



N95 CUP SHAPE

This surgical mask helps provide respiratory protection against certain airborne biological particles. It is disposable and fluid resistant to splash and spatter of blood and other infectious material.

This healthcare product is designed to help provide respiratory protection for the wearer.

Certification:

NIOSH approved N95,
TC-84A-7228.
99% BFE (Bacterial Filtration Efficiency)
GB/T 24001-2016 / ISO 14001:2015 Standard
ISO 9001:2015
ISO 45001:2018
EN 149: 2001+A1:2009



N95 FOLDABLE MASK

The foldable N95 masks are made out of 5 layers and provide over 95% bacterial filtration efficiency.

They are extremely durable with a soft and comfortable inner surface, an adjustable nose piece and secured head straps to provide a proper fit.

Certification:

NIOSH approved N95,
TC-84A-6973
95% BFE (Bacterial Filtration Efficiency)
EN 149: 2001+A1:2009
GB/T 24001-2016 / ISO 14001:2015 Standard
ISO 9001:2015
ISO 45001:2018

APPAREL



ALITMEDICAL



GENERAL ISOLATION GOWN – 25GSM

Made from lightweight polypropylene and multi-layer materials (SPP material), includes a neck tie and an elastic wrist. Suitable for isolation and basic protection from bacteria and other particles. It has strong sealed seams designed to withstand tough conditions.

Available in two sizes 115 x 137cm and 120 x 140cm, and in three colours – white, yellow and blue.

Applications:

Medical purpose / examination, healthcare, general housekeeping, laboratory, industrial purpose / PPE

Certificate:

CE approved
FDA approved



GENERAL ISOLATION GOWN – 45GSM

Made from lightweight polypropylene and multi-layer materials (SPP material), includes a neck tie and an elastic wrist. Suitable for isolation and basic protection from bacteria and other particles. It has strong sealed seams designed to withstand tough conditions.

Available in two sizes 115 x 137cm and 120 x 140cm, and in three colours – white, yellow and blue.

Applications:

Medical purpose / examination, healthcare, general housekeeping, laboratory, industrial purpose / PPE

Certificate:

CE approved
FDA approved



ISOLATION GOWN – AAMI LEVEL 2

The Isolation Gown is made from lightweight SMMS fabric and includes a neck tie and knitted cuffs.

This product is AAMI Level 2 certified and is suitable for isolation and basic protection from bacteria and other particles.

It comes in two sizes 115 x 137cm and 120 x 140cm and is available in three different colours – white, yellow and blue.

Applications:

Medical purpose / examination, healthcare, general housekeeping, laboratory, industrial purpose / PPE

Certificate:

ANSI/ AMMI PB70-2012: level 2
AATCC 42-2017
CE approved
FDA approved



ISOLATION GOWN – AAMI LEVEL 3

Our AAMI Level 3 certified isolation gowns are made from SPP+PPE material. They are waterproof and stitching seams are sealed with waterproof tapes for better protection. Suitable for isolation and basic protection from bacteria. These gowns are very comfortable featuring a neck tie, a classic tie waist and an elastic wrist. They come in one size – 120 x 140cm and are available in three different colours – white, yellow and blue.

Applications:

Medical purpose / examination, healthcare, general housekeeping, laboratory, industrial purpose / PPE

Certificate:

ANSI/ AMMI PB70-2012: level 3
AATCC 42-2017
AATCC 127-2018
CE approved
FDA approved

Not for surgical use

APPAREL



STERILE SURGICAL GOWN – AAMI LEVEL 3

Made out of SMS non-woven fabric weighing 45g/M2. SMS performance: fluid proof and anti-static.

These disposable surgical gowns are mainly used by doctors to wear during hospital surgeries, treatment of patients and home care.

All gowns have knitted cuffs and Velcro hooks and piles.

Applications: Surgical use only

Certificate:

Meets EN13795
CE approved
FDA approved

Remarks:

AMMI LEVEL 3, LATEX FREE AND
WITH EO STERILIZATION.



SHOE COVERS

Made from SPP material weighing 35gsm, they are suitable for isolation and basic protection from various particles. Soft and lightweight with great fit and performance.

Shoe caps come in one size 160 x 410mm and in blue colour.

Certificate:

CE approved
FDA approved



BOUFFANT CAPS

Our latex free bouffant cap is made of spun-bonded polypropylene and sewn with a class one flammability rating elastic band, providing great protection and comfort.

Available in three sizes 18", 21" and 24" and in blue and white colours.

Applications:

Medical purpose / examination, healthcare, general housekeeping, laboratory, industrial purpose / PPE

Certificate:

CE approved
FDA approved



MOB CAPS

Our mob caps come in single and double elastic, both 53cm in size. Made from SPP material, these caps offer a comfortable fit, secure hair protection, minimizing the risk of contamination. The design of the elastic accommodates a variety of hairstyles.

Applications:

Medical purpose / examination, healthcare, general housekeeping, laboratory, industrial purpose / PPE

Certificate:

CE approved
FDA approved

GLOVES



ALITMEDICAL

GLOVES



DISPOSABLE NITRILE EXAM GLOVES

Synthetic nitrile material is also known as thermoplastic elastomer (TPE). These gloves are made from a powder-free material and are latex-free.

There are 100 pieces per box and 10 boxes per carton.

Gloves are between 3.5–4mm thick, and come in five sizes XS, S, M, L, XL. Available in three colours – white, blue and violet.

Applications:

Medical purpose / examination, healthcare, general housekeeping, laboratory, industrial purpose / PPE

Certificate:

Meets EN455 and EN374
Meets ASTM D6319
Meets ASTM F1671
CE approved
FDA approved



SYNTHETIC VINYL EXAM GLOVES

Our synthetic vinyl exam gloves are easier to use than standard vinyl gloves and have longer wearing time than nitrile gloves. They are PU coated with no protein, limiting allergy risk (no latex proteins) and offering an increased protection.

Gloves come in five sizes – XS, S, M, L, XL – and are available in three colours – white, blue and violet.

Applications:

Medical purpose / examination, healthcare, general housekeeping, laboratory, industrial purpose

Certificate:

Meets EN455 and EN 374
Meets ASTM D5250
Meets ISO374
CE approved
FDA approved

SYRINGE & NEEDLE



ALITMEDICAL

SYRINGE & NEEDLE



STERILE SINGLE USE HYPODERMIC SYRINGE WITH / WITHOUT NEEDLE

Our single use syringe specifically designed for the curative medications.

It is manufactured with medical grade materials in compliance with International Standards that ensures its safety and effectiveness to its purpose, and embedded with best quality features for best patient comfort.

Applications:

Medical purpose / examination,
healthcare, general
housekeeping, purpose
laboratory, industrial / PPE

Certificate:

ISO.13485:2016 ISO.9001:2015
MDD-93/42/EEC CE Marking
FDA approved



STERILE SINGLE USE HYPODERMIC NEEDLE

A single use hypodermic needle manufactured using state of the art technology from Europe automation expert, KHALE.

It is produced with non-toxic and medical grade materials for its safety and effectiveness to its purpose.

The ultra-sharp needle tip is optimised for minimal penetration force thereby reducing patient discomfort

Certificate:

ISO.13485:2016 ISO.9001:2015
MDD-93/42/EEC CE Marking
FDA approved

SYRINGE & NEEDLE



STERILE SINGLE USE SYRINGE FOR INSULIN USE

Our Insulin Syringe serving as a tool in diabetic care, unique combination of ultra-fine needle with its tip optimised for minimal penetration force and variety of needle sizes.

It is manufactured on a technologically advanced machine from Europe, available in U-100 and U-40 Insulin units Certificate:

Certificate:

ISO.13485:2016 ISO.9001:2015
MDD-93/42/EEC CE Marking
FDA approved



STERILE SINGLE USE INTRAVASCULAR CATHETER / CANNULA

It is designed with specific care for infusion therapy and available with many customised variables to cater the need of healthcare professionals.

It is manufactured on a state-of-the-art technology, fully computerised and automated machines originated from Europe. It is guaranteed with its best quality with no manual operations throughout.

Certificate:

ISO.13485:2016
ISO.9001:2015 MDD-93/42/EEC
CE Marking FDA approved

SYRINGE & NEEDLE



STERILE SINGLE USE SAFETY SYRINGE WITH RE-USE PREVENTION FEATURES (RUP)

Each year, millions of healthcare workers get infected with needle stick injury (NSIs).

Needle stick injuries are among the most prevalent and preventable occupational accidents.

The ultimate solution to prevent needle stick injury is our “RUP”, the re-use prevention syringe.

It is cleverly designed for the healthcare workers, easy to use, without any additional training.

Our “RUP” is a patented product and is welcomed by many healthcare organisations.

Certificate:

ISO.13485:2016 ISO.9001:2015

MDD-93/42/EEC CE Marking

FDA approved

ALITMEDICAL



Head office: 407-8, Hing Yip Commercial Building, 284
Des Vouex Road Central, Hong Kong, HK
Tel: +852 9811 8857

UK address: 39-40 Queen's Gate, Kensington, London,
SW7 5HR, UK | Tel: +44 208 638 8232

USA address: 145 Talmadge Road, Suite 4, Edison, NJ
08817, USA | Tel: +1 732 734 0699

UAE address: 297 Bayan Business Center, Green
Community Village, Dubai, UAE | Tel: +971501355145

Email us on info@alitmedical.com

